



CERTUS Magnetic Stirrer System

Specifications:

Fluid Supply	Gyger glass syringes: 10 ml & 20ml
Mode of Use	Accessory of the CERTUS FLEX
Valve Type	Various Gyger SMLD 300GC / 300GCR
Control Type	Manual control
Motor Specs	Brushless DC
Power Supply Output	24VDC
Power Supply Input	120/220V
Frequency	50/60 Hz
Wetted material	PEEK/borosilicate glass/FFKM/FFPM
Speed	Infinitely variable from 0 -1100 RPM
Drive System	Belt driven rotary stirrer

The brand-new **Magnetic Stirrer System** for the **CERTUS FLEX** offers many advantages:

- Reusable high-quality borosilicate Gyger glass syringes
- Simple to operate
- Speed control
- Suitable for dispensing cells, beads & other particles
- Keeps particles in suspension
- Significantly reduces sedimentation and guarantees homogenous dispensing

